Outer

pigeon

Electric
Finishing Toothbrush
Cepillo de dientes electrico

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EN

WARNING Users of orthodontics & implants, please use under the guidance of a paediatric dentist. This product is suitable to be used by caregiver on child age 12 months onwards. This product is only designed to clean teeth. If your child experience any discomfort or pain during use, please stop using it and consult a paediatric dentist. This appliance is not intended for use by children under age three or persons with reduced physical or mental capabilities, unless they are given supervision by a person responsible for their safety. This is not a toy. Adult supervision always required. Do not allow children to touch the battery & screw. Must be handled by an adult. Do not adjacent to or stacked with other electrical equipment should be avoided because it could result in improper operation. Please dispose this toothbrush in accordance with your local regulations.

Caution If your child is currently suffering from tooth decay, ulcers or any other dental or gum conditions, seek your paediatric dentist's advice before using the toothbrush. Consult your paediatric dentist immediately if excessive bleeding occurs after use.

Prohibited Do not use this product while bathing. Do not place or store this product where it can fall or be pulled into a tub or sink. This product should never left unattended when turned on. Use only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer. Never operate this product if it has a damaged brush head, if it is not working properly, if it has been dropped or damaged. Do not scratch, modify, place near a heater, bend excessively, twist, pull, place heavy objects on, or bundle it. Doing so may cause failure.

Parts Name

Brush Head**
Power Button On / Off
Motor cover / Handle
Screw (to secure the Battery Cover & battery inside the Motor)
Battery Cover (inside the Motor unit)

Assemble and Disassemble / Replace the Brush Head

Disassemble
Loosen the Brush Head by twist it to the right (1). Remove it from Motor unit (2).

Assemble
Align the Brush Head to the alignment point on the Motor unit and attach the Brush Head (1) to the Motor unit. Tighten the Brush Head by twist it to the left (2) until you hear a soft 'click'.

To Insert Battery

You will need: 1x small phillips head screwdriver (not included) / 1x AAA battery (not included)

Step 1 To open, hold the Motor unit. Turn and twist Motor Cover clockwise. Pull out the Motor Cover.

Step 2 Unscrew Battery Cover

Step 3 Gently pull this part to take out the Battery Cover.

Step 4 Insert AAA+ battery

Step 5 Reassemble the Battery Cover and other part

Step 6 Push up Motor Cover. To lock, hold the Motor unit. Turn and twist Motor Cover anticlockwise.

Ready to use!

Daily Maintenance & Care

1. Please clean the product before using. Or cleaning the product every 6 month if you don't use for a long time.
2. Clean the product with clean water. Do not use corrosive or abrasive detergent. Do not use hot water over 50°C in order to avoid damaging the product. Keep out of reach of children.
3. Do not attempt to perform any servicing, maintenance or repair while the toothbrush is switched on.

Cleaning the Motor unit:

- 1. Please clean the Motor unit with a clean & soft cloth.
- 2. Do not wash the Motor unit in the water.

Cleaning the Brush Head:

- 1. After each use, rinse the Brush Head under running water (1) or Remove the Brush Head and rinse Brush Head and Handle separately (2). Then air dry.
- 2. Do not share brush.
- 3. Never put in dishwasher.

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Operation

OFF

Press 1 **LED ON** Auto-off after 1 min if inactivity.

Press 2 **CLEANING MODE 1** A slower speed mode 12,000 times/min for a gentler brushing experience but still providing a thorough clean. Spend 30 sec for each quadrant of the mouth (upper right, upper left, lower right and lower left).

Press 3 **CLEANING MODE 2** Speed mode 16,000 times/min cleaning mode for enhanced cleaning results. Spend 30 sec for each quadrant of the mouth (upper right, upper left, lower right and lower left).

Press 4 **OFF** 30 sec - "Beep" & LED Light bling once; 60 sec - "Beep Beep" & LED Light bling once; 90 sec - "Beep Beep Beep" & LED Light bling once. Pause in the brushing at 30 second intervals and then at the end of the 2 minute cycle.

Not working properly?

Not working	<input type="checkbox"/> Is it switched on? <input type="checkbox"/> Are the batteries inserted in the correct positive '+' and negative '-' direction? <input type="checkbox"/> Are the batteries out of charge?
LED Light is not working	<input type="checkbox"/> Are the batteries out of charge?
No vibration after switched on	<input type="checkbox"/> Pressing the Brush Head bristle hard on the teeth? <input type="checkbox"/> Are the batteries out of charge? <input type="checkbox"/> Is the Brush Head assembled correctly? <input type="checkbox"/> Is the bristle has worn out?

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